

BOARD OF
BUILDING AND SAFETY
COMMISSIONERS

MARSHA L. BROWN
PRESIDENT

VAN AMBATELOS
VICE-PRESIDENT

VICTOR H. CUEVAS
HELENA JUBANY
ELENORE A. WILLIAMS

CITY OF LOS ANGELES
CALIFORNIA



ANTONIO R. VILLARAIGOSA
MAYOR

DEPARTMENT OF
BUILDING AND SAFETY
201 NORTH FIGUEROA STREET
LOS ANGELES, CA 90012

ROBERT R. "BUD" OVROM
GENERAL MANAGER

RAYMOND S. CHAN, C.E., S.E.
EXECUTIVE OFFICER

Ernie Patrick
Dentalez
2500 Hwy 31 South
Bay Minette, AL 36507

Telephone: 800-383-4636

Models:

Galaxy / Advantage / Simplicity: CMU-A, -B, -C, -D, -AE, -BE, -C, -DE.

File Number: M-960090

Effective Date: 01/20/2012

Expiration Date: 01/01/2013

Manufacturer: SAME

Item: DENTAL UNITS W/WO CUSPIDORS

PLEASE REFER TO THE APPLICATION NUMBER ON ALL CORRESPONDENCE, PHONE CALLS & SAMPLES

The item processed under the above file number has been tested, examined and found to comply with applicable Los Angeles Plumbing Code (LAPC) and/or Los Angeles Mechanical Code (LAMC) and is hereby approved unless revoked for cause.

This letter or the letterhead copy may not be reproduced for use in labeling or advertising, but copies of this letter may be used within your organization and shown as evidence of approval.

CONDITIONS OF APPROVAL:

- 1. A pressure type vacuum breaker shall be installed on each water inlet connection to the delivery unit, minimum 12" above all downstream piping and flood level rim of receptor per section 94.603.4 of LAPC, 2011 Edition.**
 - 2. The drain line shall be routed without dips or sags, minimum of 1" or two pipe diameter above an approved indirect waste receptor per Section 94.801 of LAPC, 2011 Edition or shall be connected to an approved vacuum system.**
- Note: The bowl wash and cup fill both have air gaps and do not require additional backflow preventers.**

Each item shall have the name or trademark of manufacturer and model number where it will be visible for inspection.

In order to renew this approval, an application and reexamination fee must be submitted on or before the expiration date shown above. Arrangement shall be made to remove any test samples remaining in the Laboratory. Samples remaining after 21 days from the date of this notice will be shipped collect by the most convenient carrier or disposed of in accordance with applicable regulation.

Thomas Liu
Mechanical Testing Laboratory, Test Engineer

Plumbing / Mechanical Inspection
Inspection Bureau