

5400

Forest
DENTALEZ

5400 Package

A Dental Chair
Engineered for
Everyday Excellence

*Practical.
Professional.
Proven.*



5400 CHAIR

Value You Can Depend On

PRACTICAL



Clinical-First, Comfort-Focused Design

- **Sit or stand—it's your choice.** As the highest-rising chair in the industry (up to 36.5"), the 5400 delivers exceptional treatment flexibility while supporting optimal clinician comfort across procedures.
- **Industry-leading lift capacity.** Supporting up to a 525 lb patient weight (650 lb total system capacity, excluding delivery unit and light), the 5400 is a truly comprehensive chair designed to enhance stability, safety, and patient comfort.

PROFESSIONAL



Premium Construction

- **Next-generation core platform.** Building on field-tested engineering, the 5400 delivers a focused, high-performance design, engineered for durability, reliability, and long-term performance.
- **Assembled in the USA** to ensure superior quality with shorter lead times, and **exceptional value.**

PROVEN



Clinically Tried and True

- Backed by a long history of clinical performance and continuous innovation.
- For more than 100 years, DENTALEZ has delivered the ultimate dental experience—earning the trust of clinicians while enhancing patient comfort.

Features at a Glance

Essential Features for High-Volume Practices

- Fully adjustable chair positioning for clinical precision
- Ergonomic design for clinician comfort
- High-quality, durable upholstery for patient satisfaction
- Robust base and frame engineered for longevity
- Smooth, quiet operation for a calm clinical environment



Designed Around Your Practice

Proven to enhance workflow, improve patient experience, and provide consistent, reliable performance—every day.

“From a day-to-day clinical standpoint, the chair presents as solid, practical, and well-suited for high-volume, real-world dentistry. The overall design supports efficient operatory flow and ergonomics that matter during long procedures and full schedules.”

Dr. Benjamin Gates
Coeur d'Alene, ID



5400 Chair Package Includes:



5400 Patient Chair

- Fully adjustable chair positioning for clinical precision
- Ergonomic design for clinician comfort
- High-quality, durable upholstery for patient satisfaction



9077 LED Light

- Hands-free operation with automatic on/off functionality
- 3-Axis rotation to position the light exactly where needed for treatment
- 8 Energy efficient LED bulbs to create uniform light in the oral cavity

4486PI Delivery Unit

- Provides easy operation with simple integrated holders and an easy-to-read pressure gauge
- Accessible control unit for simple maintenance and less downtime



Additional Options:

Simplicity Doctor Stool - Black



Simplicity Assistant Stool - Black



Optional Footswitch - Shadow




Rear Assistant Instrumentation



DENTALEZ partners with industry leaders to drive innovation and improve outcomes for oral health professionals and their patients. We are the leader in Smart Technology.

DENTALEZ – your trusted clinical, and innovation partner, leading the way in Smart Technology in the dental industry.



Practical, Professional, & Proven. The 5400 Package delivers the essential features and performance that our clinicians rely on. Built on a foundation of premium engineering, it balances durability, patient comfort, and clinical functionality—at just the right price.

Take the Next Step: Schedule a demo today and see why **the 5400 is the chair of choice for clinicians who value quality, reliability, and patient satisfaction.**

SPECIFICATIONS

Specification	Description
Seat Height	17" - 36.5"
Upholstery	Seamless
Total Lift Capacity	525 pounds (650 patient only, less unit and light)
Lift System	Hydraulic
Headrest	Standard Dual Articulating
Armrests	Swing down
Swivel/Rotation	30 degree left/right, 60 degree total
Standard Controls	Standard touchpads on both sides of upper backrest
Control Functions	3 auto programs plus auto return and last position
Limited Warranty	5 years

DETAILS


Category	Part Number	Description
Chair	5400UC	Chair, Forest 5400, Naugahyde Upholstery
Chair	5400UL	Chair, Forest 5400, UltraLeather Upholstery
Delivery Unit	4486PI	Pivot Chair Mount Unit w/Internal Umbilical & Standard IC Control Head
Delivery Unit	4296PI	Pivot Chair Mount Unit w/Internal Umbilical & Fusion Control Head
Post	1113-654	Light, Pivot Post
Delivery Unit	1095	Assistant Instrumentation w/Internal Umbilical, Rear Pivot Chair Mount
Delivery Unit	1096	Assistant Instrumentation w/Internal Umbilical, w/Touchpad, Rear Pivot Chair Mount
Light	9077PT	Operator LED Light, Chair Mount, 5400
Optional Accessory	3680-500	Option, Footswitch, 5400, Shadow
Optional Accessory	6414-046	Option, Unit Mounted Touchpad

DENTALEZ

Connecting innovation to care

DENTALEZ is a proven leader, dedicated to providing innovative and trusted products, service, and support. We empower healthcare professionals to maximize their operations for increased patient care.

©2026 DENTALEZ, Inc. DENTALEZ is a registered trademark of DENTALEZ, Inc.
UltraLeather is a registered trademark of Ultrafabrics, Inc. Rev. 09/2026
0425 2.5M 350



Contact a Sales Rep