
E900

Star
DENTALEZ

Star E900 Electric Handpiece System

Introducing an ***entirely new*** electric handpiece system from DENTALEZ.





Smart Integration, One Comprehensive System

Track and monitor your system with RFID technology and Aeras Intel

E900

Made in Lancaster, PA, USA

Real Benefits for Every User



One Comprehensive System

- All in one unit with three factory settings for restorative and endo procedures consolidating multiple functions into a single system.
- Helps reduce chair time and increases operatory efficiency.
- Ten profile settings available where runtime hours can be viewed for each profile on the web application user interface.



Ergonomically Designed Handpieces

- Helps reduce hand and wrist fatigue with comfortable and balanced weight distribution.
- The Star E900 handpiece attachments are shorter and 15% lighter* in weight than competitive handpieces.
- Designed to reduce tubing drag on wrist.


*Competitive analysis based on an internal study, data on file.



Actively Monitors Operating Temperature

- A sensor in the handpiece** detects the high temperature, which then notifies the user with a screen notification and shuts off the motor.
- Notifies the user if a motor replacement is recommended.

**Available on Star E900 1:5 attachment.



Protect Your Investment with Smart Technology

- Manage and track your assets^ with built in RFID in the handpiece, motor and tubing.
- Know where your assets are in real time.

^With compatible reader system.


Made in Lancaster, PA, USA

Use the Aeras web app to track the handpiece hours and performance. You can choose, customize and save the user profiles on the touchscreen on the Star E900 Electric System.






Compatible with short coupling (nose lengths)



Prep Mode Screen



Endo Mode Screen

Made in Lancaster, PA, USA

The Star E900 Electric Handpiece System by DENTALEZ:
Work Longer. Work Safer. Work Smarter.
Engineered for excellence inspired by you.

Part Number	Description
266340	Star E900 Electric System, Sterling Tubing
266341	Star E900 Electric System, Black Tubing
266345	Star E900 1:1 Straight Handpiece
266343	Star E900 1:1 Auto Latch, Fiber Optic, Handpiece
266342	Star E900 1:5 Auto Chuck, Fiber Optic, Handpiece
266490	Star E900 1:5 Auto Chuck, Non Fiber Optic, Handpiece



Compatible with competitive systems — attachments plug into your current set-up.
Temperature sensing applies only when the Star E900 1:5 attachment is connected to the Star E900 Electric System.



Star E900 1:5
Auto Chuck, Fiber Optic, Handpiece



Star E900 1:1
Auto Latch, Fiber Optic, Handpiece



Star E900 1:1
Straight Handpiece

DENTALEZ

Connecting innovation to care

©2026 DENTALEZ, Inc. All rights reserved. DENTALEZ, and Star are registered trademarks of DENTALEZ INC. AERAS, Inc. is a separate corporation from DENTALEZ, Inc.

E900

Made in Lancaster, PA, USA